




Transonic®
Surgical Flowmeters

HT313/HT323 (-CS) Specifications


Intraoperative Blood Flow Measurements

- ◆ Cardiothoracic
- ◆ Vascular
- ◆ Transplant
- ◆ Neurosurgery

“Push and Play” ease of operation with the essential features for measuring and documenting blood flow during surgery



The HT313 & Portable Trolley
A 2-meter extension cable allows positioning at the most convenient location in the OR.





Print FlowSound™ Invert Signal Quality Indicator Mean Flow

Easy to Operate Plug in flowprobe & measure

Large Display Mean flow in ml or L/min can be read across room.

Documentation A flow recording printout for your patient's chart.

FlowSound™
 Listen to FlowSound while focus remains on the surgical field. Pitch equals volume flow.

 **Signal Quality Indicator** Indicates quality of ultrasonic contact.

Invert Changes a negative reading to a positive reading so that the flowprobe can be applied to either side of the vessel.

Central Monitor Interface
Optional interface to central patient monitor station for continuous real-time display.



HT313/HT323 (-CS) Specifications

PHYSICAL SPECIFICATIONS

single/dual channel flowmeters with built-in chart recorder;
HT313(-CS): 10 lb. (4.5 Kg); HT323(-CS): 11 lb. (5 Kg);
HT313(-CS)/HT323(-CS): 11" (27.8 cm) wide x 4" (9.12 cm)
high x 14" (35.4 cm) deep; attached two meter flowprobe-to-
flowmeter extension cable per channel, separate grounded
medical grade power line-cord is supplied with the flowmeter.

INPUT POWER

Universal Power Supply; 50-60 Hz; 100-240 V \pm 10%
(automatic voltage adjustment)

FUSES

Both power entry conductors are protected by fuses (5 x 20
mm fuses) 1.5A fast blo, mfg: bussman #GMA1.5, 250 VAC)

ELECTRICAL ISOLATION

Flowmeter cabinet is grounded; line to ground leakage
current: less than 50 μ A @100-120 V line, less than 100 μ A @
220-240 V line. All electronic components and cabling of the
probe extension cable and probe are fully isolated to meet IEC
"cardiac floating" specifications. Defibrillator protection: probes
can remain attached to the patient during cardiac defibrillation
for instantaneous report on restoration of flow.

AUTOMATIC METER ADJUSTMENTS

- Ultrasound frequency and insonification parameters
- Probe size and corresponding flow output ranges
- Volume flow calibration of connected flowprobe
- Flow sampling rate optimized for local acoustic conditions

DIGITAL IDENTIFICATION

Probe identification and calibration parameters are
programmed on an EEPROM inside the probe connector.

FLOW MONITOR ANALOG OUTPUTS

Uniscale Output: 10 Hz low pass filtered; automatic gain
adjustment for fixed scaling of multiple probe sizes. Accepts
mini-phono connector. HT313: 1 uniscale output; HT323: 2
uniscale outputs.

AUXILIARY (-CS) OPTION FLOWMETER INPUTS

Two BNC connectors are provided on the rear panel as
auxiliary input ports (i.e., pressure or EKG monitors).
Each accepts a signal range of \pm 5 Volts.
This option includes the FlowTrace software package.

FLAWSOUND™

Audio representation of flow sensed by the probe; one octave
increase in sound pitch = 4x increase in volume flow.

ULTRASOUND SIGNAL QUALITY CHECK

Front panel ultrasound "container-style" level indicator provides
visual feedback on ultrasound signal strength; automatically
disables flow measurements when ultrasound signal conditions
are insufficient. Front panel 6-digit display reports error
condition.

CHART RECORDER

Prints a permanent record of pulsatile flow wave along with
date stamp and values for mean, peak and minimum flow.
Front panel push button engages the recorder; rear panel
switch selects from 3 printing modes: "Slow" prints 20 cm flow
record at 20 mm/sec; "Fast" prints 20 cm flow record at 40
mm/sec; "RUN/STOP" prints a continuous strip at 20 mm/sec.

COMPATIBLE ULTRASONIC TRANSDUCERS

These meters accept a range of HQ-Series probes for
intraoperative surgical use. See probe specification sheets or
your local representative for models, sizes and recommended
uses.

ULTRASONIC FREQUENCY/PARAMETERS

The ultrasound output level of the flowsensors is factory-set
and does not incorporate any interactive system features.
These settings are made using "ALARA" principles (As Low As
Reasonably Achievable), and are orders of magnitude below the
FDA "preamendment levels," the USA insonification safety
limits. All are probe size dependent.

Transducer excitation: Burst of 10 to 24 waves
Transducer excitation frequency: 1.2 MHz to 4.8 MHz
Transducer excitation rate: 1.2 KHz to 14 KHz
automatically reduced further when probe operates in adverse
acoustic conditions that would degrade measurement
accuracy.

USB-PORT

Optically isolated output connector meeting IEC0601 patient
isolation. Permits connection to a laptop or other data monitor
computer. The \pm 5V analog flow signals (and pressure, EKG
Signals for -CS option meters) are digitized at 100 sampling Hz
rate with 12-bit resolution. Data is supplied at 19200 Baud for
use with FLOWTRACE software

REGULATORY COMPLIANCE

Transonic Systems flowmeters and sensors comply with USA
standards for medical and dental equipment (UL2601-1), and
with European standards for medical and ultrasonic apparatus
(IEC0601-1). These products are CE marked per
93/42EEC Annex II. Transonic Systems is an
ISO9001/EN46001-certified facility.



Flowprobe Selection Guide

PERIVASCULAR FLOWPROBE SERIES & AVAILABLE SIZES		
SUFFIX	DESCRIPTION	SIZES (mm)
-FMC	Coronary	1.5, 2, 3, 4
-FMV	Vascular	1.5, 2, 3, 4, 6, 8, 10, 12, 14
-FME	Carotid (L-reflector for carotid endarterectomy)	4, 6, 8, 10
-FTV	OptiMax® (hands-free, J-reflector)	4, 6, 8, 10, 12
-FTE	OptiMax® (hands-free, L-reflector carotid endarterectomy)	4, 6, 8, 10, 12
-FMU	Microvascular (handle)	0.7, 1, 1.5, 2, 3
-AU	Cardiac Output COnfidence Flowprobe®	8, 10, 12, 14, 16, 20, 24, 28, 32, 36
-MB	Intracranial Charbel Micro-Flowprobe® (long bayonet handle)	1.5, 2, 3
-MB-S	Extracranial EC-IC Bypass: Micro-Flowprobe (short bayonet handle)	3, 4, 6
-FSB	Basic (no handle, L-reflector, sliding cover)	1.5, 2, 3, 4, 6, 8, 10, 12, 14

Recommended Sizes and/or Flowprobe Series for Specific Vessels or Applications

CARDIAC SURGERY		
CABG: ON OR OFF PUMP	Probe Size (mm)	Probe Series
Arterial conduits	1.5, 2, 3, 4	-FMC
Saphenous vein	2, 3, 4	-FMC
CARDIAC OUTPUT		
Ascending aorta	28, 32, 36	-AU
Pulmonary artery	24, 28, 32	-AU
Pediatric heart	8, 10, 12, 14, 16, 20	-AU

TRANSPLANT SURGERY		
LIVER	Probe Size (mm)	Probe Series
Hepatic artery	4, 6, 8	-FMV -AU
Portal vein	10, 12, 14	-FMV -AU
KIDNEY		
Renal artery	4, 6	-FMV -FSB
Renal vein	10	-FMV -FSB
External iliac artery	6, 8	-FMV -FSB
Hypogastric artery	4, 6	-FMV -FSB
PANCREAS		
Common iliac artery	8	-FMV -FSB

CEREBROVASCULAR SURGERY		
ANEURYSM CLIPPING	Probe Size (mm)	Probe Series
Cerebral arteries	1.5, 2, 3	-MB -MR
EC-IC BYPASS		
Extracranial	3, 4, 6	-MB-S MR-S
Intracranial	1.5, 2, 3	-MB -MR
AVM, TUMOR RESECTION, DURAL FISTULA		

VASCULAR SURGERY				
CAROTID ENDARTERECTOMY	Probe Size (mm)	Probe Series		
Common carotid artery	8, 10	-FTE	-FME	-FSB
External carotid artery	6	-FTE	-FME	-FSB
Internal carotid artery	6	-FTE	-FME	-FSB
AV FISTULAS & GRAFTS				
Radial artery	2, 3	-FMV		-FSB
Brachial artery	3, 4, 6	-FMV	-FTV	-FSB
Graft venous outflow	4, 6	-FMV	-FTV	-FSB
ABDOMINAL				
Renal bypass	4, 6	-FMV	-FTV	-FSB
Aortoiliac shunt				
aorta	16, 20			-AU
common iliac	10, 12	-FMV	-FTV	-AU -FSB
Portocaval shunt	10, 12, 14	-FMV	-FTV	-AU -FSB
Splenorenal shunt	10, 12, 14	-FMV	-FTV	-AU -FSB
LOWER EXTREMITY BYPASS				
Profunda femoris	8	-FMV	-FTV	-AU -FSB
Common femoral	8, 10	-FMV	-FTV	-AU -FSB
Popliteal	4, 6	-FMV	-FTV	-FSB
Tibial	3, 4	-FMV	-FTV	-FSB

MICROVASCULAR SURGERY		
REATTACHMENTS/FLAPS	Probe Size (mm)	Probe Series
Microvessels in hand, wrist	0.7, 1, 1.5, 2, 3	-FMU

