

TRI-FLOW™

FEATURES & BENEFITS

- Ideal for extracorporeal therapies: dialysis, hemofiltration, apheresis and plasmapheresis.
- Third lumen gives greater flexibility for fluid administration, vasoactive drug therapy, nutrition and blood sampling.
- Integral soft tip reduces occurrence of stenosis and erosion.

Catalog #	Catheter Only	
XTP3114CT	11.5F x 12CM, STRAIGHT TRI-FLOW™ CATHETER	10/Box
XTP3116CT	11.5F x 15CM, STRAIGHT TRI-FLOW™ CATHETER	10/Box
XTP3118CT	11.5F x 20CM, STRAIGHT TRI-FLOW™ CATHETER	10/Box
XTP3119CT	11.5F x 24CM, STRAIGHT TRI-FLOW™ CATHETER	10/Box
XTP3114CTPC	11.5F x 12CM, PRE-CURVED TRI-FLOW™ CATHETER	10/Box
XTP3116CTPC	11.5F x 15CM, PRE-CURVED TRI-FLOW™ CATHETER	10/Box
XTP3118CTPC	11.5F x 20CM, PRE-CURVED TRI-FLOW™ CATHETER	10/Box

CATHETER ONLY CONTAINS:

(1) Catheter (1) Dilator (3) Injection Caps (1) Patient I.D. Card

Catalog #	Basic Set	
XTP3114MT	11.5F x 12CM, STRAIGHT TRI-FLOW™ CATHETER	5/Box
XTP3116MT	11.5F x 15CM, STRAIGHT TRI-FLOW™ CATHETER	5/Box
XTP3118MT	11.5F x 20CM, STRAIGHT TRI-FLOW™ CATHETER	5/Box
XTP3119MT	11.5F x 24 CM, STRAIGHT TRI-FLOW™ CATHETER	5/Box
XTP3114MTPC	11.5F x 12CM, PRE-CURVED TRI-FLOW™ CATHETER	5/Box
XTP3116MTPC	11.5F x 15CM, PRE-CURVED TRI-FLOW™ CATHETER	5/Box
XTP3118MTPC	11.5F x 20CM, PRE-CURVED TRI-FLOW™ CATHETER	5/Box

BASIC SET CONTAINS:

(1) Catheter (1) Introducer Needle (1) Guidewire (1) Dilator (1) Adhesive Wound Dressing (3) Injection Caps (1) Patient I.D. Card

Catalog #	Full Tray	
MCDLT3-3500	11.5F x 12CM STRAIGHT TRI-FLOW™ CATHETER	5/Box
MCDLT3-4000	11.5F x 15CM STRAIGHT TRI-FLOW™ CATHETER	5/Box
MCDLT3-4400	11.5F x 20CM STRAIGHT TRI-FLOW™ CATHETER	5/Box
MCDLT3-4424	11.5F x 24CM STRAIGHT TRI-FLOW™ CATHETER	5/Box
MCDLT3-3500PC	11.5F x 12CM PRE-CURVED TRI-FLOW™ CATHETER	5/Box
MCDLT3-4000PC	11.5F x 15CM PRE-CURVED TRI-FLOW™ CATHETER	5/Box
MCDLT3-4400PC	11.5F x 20CM PRE-CURVED TRI-FLOW™ CATHETER	5/Box

FULL TRAY CONTAINS:

(1) Catheter (1) Introducer Needle (1) Guidewire (1) Dilator (1) Safety Scalpel (1) Adhesive Wound Dressing (3) Injection Caps (8) 4"x4" Gauze (2) Chlorhexidine (2) 5cc Syringe (3) Pre-Filled Saline Syringes (1) 2-0 Silk Suture (3) Safety Hypodermic Needles (1) Hemostat (1) Fenestrated Drape (1) 5cc Lidocaine (1) Filter Straw (1) Surgical Gloves (1) Needle Stick Pad (1) Patient I.D. Card



TekMed

1300 720 727

12 Bellevue Crescent
Preston Victoria 3072 Australia

T (03) 9487 9999 F (03) 9484 8792

info@tekmed.com.au • www.tekmed.com.au

USA • ISO 13485:2003

Excerpt From:
PN2114 Rev. H 4/10