

DUO-FLOW®

FEATURES & BENEFITS

- Thermosensitive polyurethane material softens after insertion for increased patient safety.
- Raulerson IJ Catheter has a unique pre-formed curve that overcomes patient discomfort and significantly decreases the risk for infections, compared with straight catheters.
- Coaxial catheter design allows for improved flow rates and outcomes, when compared with standard double “D” dual lumen catheters.

Catalog #	Basic Set	
	Straight Catheter	
XTP94MT	9F x 12CM STRAIGHT DUO-FLOW® CATHETER	5/Box
XTP96MT	9F x 15CM STRAIGHT DUO-FLOW® CATHETER	5/Box
XTP98MT	9F x 20CM STRAIGHT DUO-FLOW® CATHETER	5/Box
	Pre-Curved Catheter	
XTP94IJS	9F x 12CM RAULERSON IJ DUO-FLOW® CATHETER	5/Box
XTP96IJS	9F x 15CM RAULERSON IJ DUO-FLOW® CATHETER	5/Box
XTP98IJS	9F x 20CM RAULERSON IJ DUO-FLOW® CATHETER	5/Box

BASIC SET CONTAINS:

- (1) Catheter (1) Dilator (1) Guidewire
- (1) Adhesive Wound Dressing (1) Introducer Needle
- (2) Injection Caps

Catalog #	Basic Set	
	Straight Catheter	
XTP114MTE	11.5F x 12CM STRAIGHT DUO-FLOW® CATHETER	10/Box
XTP116MTE	11.5F x 15CM STRAIGHT DUO-FLOW® CATHETER	10/Box
XTP118MTE	11.5F x 20CM STRAIGHT DUO-FLOW® CATHETER	10/Box
XTP119MTE	11.5F x 24CM STRAIGHT DUO-FLOW® CATHETER	10/Box
	Raulerson IJ Catheter	
XTP114IJS	11.5F x 12CM RAULERSON IJ DUO-FLOW® CATHETER	10/Box
XTP116IJS	11.5F x 15CM RAULERSON IJ DUO-FLOW® CATHETER	10/Box
XTP118IJS	11.5F x 20CM RAULERSON IJ DUO-FLOW® CATHETER	10/Box
	Pre-Curved Catheter	
XTP114IJS-2	11.5F x 12CM PRE-CURVED DUO-FLOW® CATHETER	5/Box
XTP116IJS-2	11.5F x 15CM PRE-CURVED DUO-FLOW® CATHETER	5/Box
XTP118IJS-2	11.5F x 20CM PRE-CURVED DUO-FLOW® CATHETER	5/Box

BASIC SET CONTAINS:

- (1) Catheter (1) Dilator (1) Guidewire (1) Introducer Needle
- (2) Injection Caps

US and Foreign Patents Pending.

SOFT-LINE®

Catalog #	Basic Set	
	Straight Catheter	
T73M	7F x 7CM STRAIGHT SOFT-LINE® CATHETER SET	5/Box
T74M	7F x 10CM STRAIGHT SOFT-LINE® CATHETER SET	5/Box
T94M	9F x 12CM STRAIGHT SOFT-LINE® CATHETER SET	5/Box
T96M	9F x 15CM STRAIGHT SOFT-LINE® CATHETER SET	5/Box
T98M	9F x 20CM STRAIGHT SOFT-LINE® CATHETER SET	5/Box
	Pre-Curved Catheter	
T114IJS-2	11.5F x 12CM PRE-CURVED SOFT-LINE® CATHETER SET	5/Box
T116IJS-2	11.5F x 15CM PRE-CURVED SOFT-LINE® CATHETER SET	5/Box
T118IJS-2	11.5F x 20CM PRE-CURVED SOFT-LINE® CATHETER SET	5/Box

BASIC SET CONTAINS:

- (1) Catheter (1) Dilator (1) Guidewire
- (1) Adhesive Wound Dressing
- (1) Introducer Needle (2) Injection Caps

Catalog #	Basic Set	
T114ME	11.5F x 12CM STRAIGHT SOFT-LINE® CATHETER SET	10/Box
T116ME	11.5F x 15CM STRAIGHT SOFT-LINE® CATHETER SET	10/Box
T118ME	11.5F x 20CM STRAIGHT SOFT-LINE® CATHETER SET	10/Box

BASIC SET CONTAINS:

- (1) Catheter (1) Dilator (1) Guidewire (1) Introducer Needle
- (2) Injection Caps

