



## HEMO-CATH® LT

### FEATURES & BENEFITS

- Silicone material provides biocompatibility and kink resistance.
- Ideal for extracorporeal therapies: dialysis, hemofiltration, apheresis and plasmapheresis.
- Radiopaque silicone material.

Catalog #	Straight Catheter Basic Set	
SL18P	8F x 18CM SILICONE HEMO-CATH® BASIC SET	5/Box
SL24P	8F x 24CM SILICONE HEMO-CATH® BASIC SET	5/Box

**BASIC SET CONTAINS:**

- (1) Catheter (1) Tunneler (1) Peelable Introducer (1) Guidewire
- (1) 10cc Syringe (1) Adhesive Wound Dressing (1) Catheter Clip
- (1) Introducer Needle (2) Injection Caps

Catalog #	Straight Catheter Basic Set	
MC101241	12.5F x 15CM SILICONE DOUBLE LUMEN	5/Box
MC101242	12.5F x 18CM SILICONE DOUBLE LUMEN	5/Box
MC101243	12.5F x 24CM SILICONE DOUBLE LUMEN	5/Box

**BASIC SET CONTAINS:**

- (1) Catheter (2) Injection Caps (1) Scalpel (1) Valved Peelable Introducer
- (1) Introducer Needle (1) Guidewire (1) Adhesive Wound Dressing
- (1) Catheter Clip (1) Trocar

Catalog #	Straight Catheter Basic Set	
SL28E.	12.5F x 28CM SILICONE DOUBLE LUMEN	5/Box
SL32E.	12.5F x 32CM SILICONE DOUBLE LUMEN	5/Box
	<b>Pre-Curved Catheter Basic Set</b>	
SL28PCE.	12.5F x 28CM PC SILICONE DOUBLE LUMEN	5/Box
SL32PCE.	12.5F x 32CM PC SILICONE DOUBLE LUMEN	5/Box

**BASIC SET CONTAINS:**

- (1) Catheter (2) Injection Caps (1) Peelable Introducer (1) Introducer Needle
- (1) Guidewire (1) Trocar (1) Catheter Clip (Straight Catheter Only)