

HEMO-FLOW®

FEATURES & BENEFITS


- High flow rates in excess of 400ml/min.
- Radiopaque polyurethane material that is both alcohol and iodine compatible.
- Available in straight and pre-curved designs to accommodate patient need and comfort.
- Economical.

Catalog #	Straight	
HFS24E.	14.5F x 24CM HEMO-FLOW® SET	5/Box
HFS28E.	14.5F x 28CM HEMO-FLOW® SET	5/Box
HFS32E.	14.5F x 32CM HEMO-FLOW® SET	5/Box
HFS36E.	14.5F x 36CM HEMO-FLOW® SET	5/Box
HFS40E.	14.5F x 40CM HEMO-FLOW® SET	5/Box

Catalog #	Pre-Curved	
HFS24PCE.	14.5F x 24CM PC HEMO-FLOW® SET	5/Box
HFS28PCE.	14.5F x 28CM PC HEMO-FLOW® SET	5/Box
HFS32PCE.	14.5F x 32CM PC HEMO-FLOW® SET	5/Box
HFS36PCE.	14.5F x 36CM PC HEMO-FLOW® SET	5/Box

BASIC SET CONTAINS:

- (1) Catheter (1) Introducer Needle (1) Guidewire (2) Dilators (1) Tunneler
 (1) Valved Peelable Introducer (2) Injection Caps



TekMed

T (03) 9487 9999

12 Bellevue Crescent
 Preston Victoria 3072 Australia
 F (03) 9484 8792

info@tekmed.com.au • www.tekmed.com.au