



**Silicone**  
Material

**Straight,  
Coiled**  
Tip Designs


**10.5F, 15F**  
French Sizes

**Single,  
Double**  
Cuffs



**COILED TIP**

Provides a gentler, less traumatic, showering effect of dialysate on the peritoneum and prevents catheter migration.



**SINGLE OR DOUBLE CUFF**

Ensures device securement, additional protection from bacterial contamination and prevents leaks.



**ONE-PIECE FEMALE**

**LUER ADAPTER**  
Makes connections simple and manageable.



CATHETER ONLY - STRAIGHT		
MDC31S	31CM I-SERIES, STRAIGHT WITH DOUBLE CUFF	5/BOX
MDC32S	32.5CM I-SERIES, STRAIGHT WITH DOUBLE CUFF	5/BOX
MDC37S	37CM I-SERIES, STRAIGHT WITH DOUBLE CUFF	5/BOX
MDC42S	42CM I-SERIES, STRAIGHT WITH DOUBLE CUFF	5/BOX
MDC47S	47CM I-SERIES, STRAIGHT WITH DOUBLE CUFF	5/BOX
MSC31S	31CM I-SERIES, STRAIGHT WITH SUBCUTANEOUS SINGLE CUFF	5/BOX
MSC37S	37CM I-SERIES, STRAIGHT WITH SUBCUTANEOUS SINGLE CUFF	5/BOX
MSC42S	42CM I-SERIES, STRAIGHT WITH SUBCUTANEOUS SINGLE CUFF	5/BOX
MSC46S	46CM I-SERIES, STRAIGHT WITH SUBCUTANEOUS SINGLE CUFF	5/BOX
MPP41S	41CM I-SERIES, STRAIGHT WITH PREPERITONEAL SINGLE CUFF	5/BOX
MPP46S	46CM I-SERIES, STRAIGHT WITH PREPERITONEAL SINGLE CUFF	5/BOX

BASIC SET - STRAIGHT		
MPD-242	42CM I-SERIES, STRAIGHT WITH DOUBLE CUFF	5/BOX
MPD-146S	46CM I-SERIES, STRAIGHT WITH SUBCUTANEOUS SINGLE CUFF	5/BOX
MPD-141	41CM I-SERIES, STRAIGHT WITH PREPERITONEAL SINGLE CUFF	5/BOX
MPD-146	46CM I-SERIES, STRAIGHT WITH PREPERITONEAL SINGLE CUFF	5/BOX

10.5F BASIC SET - STRAIGHT		
MPD-137S	10.5F X 37CM STRAIGHT PD WITH SINGLE CUFF	5/BOX

CATHETER ONLY - COILED		
MDC57C	57.5CM I-SERIES, COILED WITH DOUBLE CUFF	5/BOX
MDC63C	63CM I-SERIES, COILED WITH DOUBLE CUFF	5/BOX
MPP40C	40.25CM I-SERIES, COILED WITH PREPERITONEAL SINGLE CUFF	5/BOX
MPP57C	57.5CM I-SERIES, COILED WITH PREPERITONEAL SINGLE CUFF	5/BOX
MPP60C	60CM I-SERIES, COILED WITH PREPERITONEAL SINGLE CUFF	5/BOX

BASIC SET - COILED		
MPD-257	57.5CM I-SERIES, COILED WITH DOUBLE CUFF	5/BOX
MPD-263	63CM I-SERIES, COILED WITH DOUBLE CUFF	5/BOX
MPD-160	60CM I-SERIES, COILED WITH PREPERITONEAL SINGLE CUFF	5/BOX

ACCESSORIES		
MP217	17F VASCU-SHEATH INTRODUCER SET	5/BOX
MC50	MALE LUER LOCK END CAP	100/BOX
MC51	BETA CLAMP	100/BOX
MBCA	BARBED LUER LOCK CATHETER ADAPTOR	20/BOX
MCS	CATHETER STYLET	5/BOX
MTS	TUNNELING STYLET	10/BOX
MTC	TITANIUM CONNECTOR	5/BOX

**CATHETER ONLY CONTAINS:**

- (1) Catheter
- (1) Female Luer
- (1) Clamp
- (1) Injection Cap
- (1) PD Catheter IFU

**BASIC SET CONTAINS:**

- (1) Catheter
- (1) Female Luer
- (1) Clamp
- (1) Introducer Needle
- (1) Guidewire
- (1) Scalpel
- (1) Peelable Sheath Introducer
- (1) Tunneler
- (1) 10cc Syringe
- (1) Injection Cap
- (6) 4" x 4" Gauze
- (1) PD Catheter IFU

**10.5F BASIC SET CONTAINS:**

- (1) Catheter
- (1) Clamp
- (1) Introducer Needle
- (1) Guidewire
- (1) Scalpel
- (1) Peelable Sheath Introducer
- (1) 10cc Syringe
- (1) Injection Cap
- (6) 4" x 4" Gauze
- (1) Extension
- (1) PD Catheter IFU