



Silicone Material	Straight, Coiled Tip Designs	15F French Size	Double Cuff
-----------------------------	--	---------------------------	-----------------------

COILED TIP

Provides a gentler, less traumatic, showering effect of dialysate on the peritoneum and prevents catheter migration.




DOUBLE CUFF

Ensures optimal outcomes by allowing one cuff to be placed in the rectus muscle and the other in the subcutaneous tissue.

PREFORMED BEND

The lateral or downward direction minimizes the risk of exit-site infection.



CATHETER ONLY - STRAIGHT		
MC200301	STRAIGHT, 43CM LEFT V-NECK WITH DOUBLE CUFF	5/BOX
MC200302	STRAIGHT, 43CM RIGHT V-NECK WITH DOUBLE CUFF	5/BOX
BASIC SET - STRAIGHT		
MC200101	STRAIGHT, 43CM LEFT V-NECK WITH DOUBLE CUFF	5/BOX
MC200102	STRAIGHT, 43CM RIGHT V-NECK WITH DOUBLE CUFF	5/BOX
CATHETER ONLY - COILED		
MC200401	COILED, 63CM LEFT V-NECK WITH DOUBLE CUFF	5/BOX
MC200402	COILED, 63CM RIGHT V-NECK WITH DOUBLE CUFF	5/BOX
BASIC SET - COILED		
MC200201	COILED, 63CM LEFT V-NECK WITH DOUBLE CUFF	5/BOX
MC200202	COILED, 63CM RIGHT V-NECK WITH DOUBLE CUFF	5/BOX
ACCESSORIES		
MP217	17F VASCU-SHEATH INTRODUCER SET	5/BOX
MC50	MALE LUER LOCK END CAP	100/BOX
MC51	BETA CLAMP	100/BOX
MBCA	BARBED LUER LOCK CATHETER ADAPTOR	20/BOX
MCS	CATHETER STYLET	5/BOX
MTS	TUNNELING STYLET	10/BOX
MTC	TITANIUM CONNECTOR	5/BOX

CATHETER ONLY CONTAINS:

- (1) Catheter
- (1) Female Luer
- (1) Clamp
- (1) Injection Cap
- (1) PD Catheter IFU

BASIC SET CONTAINS:

- (1) Catheter
- (1) Female Luer
- (1) Clamp
- (1) Introducer Needle
- (1) Guidewire
- (1) Scalpel
- (1) Peelable Sheath Introducer
- (1) Tunneler
- (1) 10cc Syringe
- (1) Injection Cap
- (6) 4" x 4" Gauze
- (1) PD Catheter IFU

