



Silicone	Coiled	15F
Material	Tip Design	French Size



COILED TIP

Provides a gentler, less traumatic, showering effect of dialysate on the peritoneum and prevents catheter migration.



LARGE SIDE HOLES

Increase flow and reduce incidence of omental occlusion.

TITANIUM BARBED CONNECTOR

Securely joins the two catheter end segments for upper abdominal or Presternal placement.



CATHETER ONLY		
MC20XC	15F X 64CM X-SERIES, COILED	5/BOX

BASIC SET		
MC20XS	15F X 64CM X-SERIES, COILED	5/BOX

ACCESSORIES		
MP217	17F VASCU-SHEATH INTRODUCER SET	5/BOX
TRX	X-SERIES TUNNELING ROD	5/BOX
MTC	TITANIUM CONNECTOR	5/BOX
FSX	X-SERIES FALLER STYLET	5/BOX
40501	PREOPERATIVE MARKING STENCIL PACK	1/BOX

CATHETER ONLY CONTAINS:

- (1) Catheter
- (1) Titanium Barbed Connector
- (1) Female Luer
- (1) Clamp
- (1) Injection Cap
- (1) PD Catheter IFU

BASIC SET CONTAINS:

- (1) Catheter
- (1) Titanium Barbed Connector
- (1) Female Luer
- (1) Clamp
- (1) Injection Cap
- (1) Introducer Needle
- (1) Guidewire
- (1) Scalpel
- (1) Peelable Sheath Introducer
- (1) Tunneler
- (1) 10cc Syringe
- (6) 4" x 4" Gauze
- (1) Injection Cap
- (1) PD Catheter IFU

ACCESSORIES		
MP217	17F VASCU-SHEATH INTRODUCER SET	5/BOX
MC50	MALE LUER LOCK END CAP	100/BOX
MC51	BETA CLAMP	100/BOX
MBCA	BARBED LUER LOCK CATHETER ADAPTOR	20/BOX
MCS	CATHETER STYLET	5/BOX
MTS	TUNNELING STYLET	10/BOX
MTC	TITANIUM CONNECTOR	5/BOX

