



DUO-SPLIT[®]

ORDERING INFORMATION

DUO-SPLIT[®] SHORT TERM HEMODIALYSIS CATHETER

STRAIGHT	CATHETER ONLY (10 / box)	BASIC SET (5 / box)	FULL TRAY (5 / box)
13FR x 12cm	DSP134C	DSP134S	DSP134T
13FR x 15cm	DSP136C	DSP136S	DSP136T
13FR x 20cm	DSP138C	DSP138S	DSP138T
13FR x 24cm	DSP139C	DSP139S	DSP139T
13FR x 30cm	DSP130C	DSP130S	-----
PRE-CURVED			
13FR x 12cm	DSP134PC	DSP134PS	DSP134PT
13FR x 15cm	DSP136PC	DSP136PS	DSP136PT
13FR x 20cm	DSP138PC	DSP138PS	DSP138PT
13FR x 24cm	DSP139PC	DSP139PS	DSP139PT
RAULERSON IJ			
13FR x 12cm	DSP134IJC	DSP134IJS	DSP134IJT
13FR x 15cm	DSP136IJC	DSP136IJS	DSP136IJT
13FR x 20cm	DSP138IJC	DSP138IJS	DSP138IJT
13FR x 24cm	DSP139IJC	DSP139IJS	DSP139IJT

Catheter Only Components

- (1) Catheter w/Dual Stylets
- (2) Injection Caps
- (1/2) Vessel Dilator(s)

Basic Set Components

- (1) Catheter w/Dual Stylets
- (2) Injection Caps
- (1) Introducer Needle
- (1) Guidewire
- (2) Vessel Dilators
- (1) Adhesive Dressing

Full Tray Components

- (1) Catheter w/Dual Stylets
- (2) Injection Caps
- (1) Safety Introducer Needle
- (3) Safety Hypodermic Needles
- (1) Guidewire
- (2) Vessel Dilators
- (1) Adhesive Dressing
- (1) Safety Scalpel
- (1) Fenestrated Drape
- (1) Hemostat
- (2) CSR Wraps
- (4) 4" x 4" Gauze Pads
- (1) Chloraprep
- (1) 5cc Syringe
- (1) 10cc Syringe
- (1) 5cc Ample Lidocaine
- (1) Set, Surgical Gloves
- (1) 2-0 Silk Sutures

Short Term
Hemodialysis
Split Tip Catheter



Innovative Designs Advancing Care Improving Lives

DUO-SPLIT®

Short Term Split Tip Dialysis Catheter

Medcomp® has once again advanced catheter technology by introducing the **Duo-Split®** Short Term Dialysis Catheter. The unique Split Tip technology with open distal tip design provides increased flow and patency rates. By incorporating a dual “Fast Track Stylet” system, healthcare professionals can follow standard insertion techniques for a fast and less invasive insertion. Finally, clinicians and patients can receive the benefits, proven from the Split Cath® long term dialysis catheter, in the short term **Duo-Split®** dialysis catheter.

SPLIT TIP

- Open tip design allows for greater flow and improved patency when compared to standard tapered tip temporary catheters.
- The first split tip short term catheter.

- Dual stylets allow standard insertion technique.

Up to 30cm catheter length

CATHETER OPTIONS

- **Lengths**
12cm
15cm
20cm
24cm
30cm

- **Configurations**
Straight
Pre-Curved
Raulerson IJ

- **Kits**
Catheter Only
Basic Set
Full Tray

Raulerson IJ

Straight

Pre-Curved

